



# **Product guide**

# Welcome

Accora is home to a team bursting with innovation and enthusiasm. We operate in the healthcare marketplace, creating and realising products for clinicians and care providers.

Our customers tell us they want products that are clinically effective at a price they can afford. By the way, we also know you want great customer service and products that go on and on and on.

## Contents

Our team	4
Configura® chairs	9
FloorBed®	27
Junior Beds	31
CommunityBed	37





Nicholas Drake Product Manager FloorBed® 01223 206100 nicholas.drake@accora.care

**Matt Hamilton** 

Account Manager 07584 857569

matt.hamilton@accora.care



Richard Smith Sales Director 01223 206100 richard.smith@accora.care

Danny Norman

Account Manager

07341 074834

danny.norman@accora.care



Daniel Smith Product Manager Configura\* 01223 206100 daniel.smith@accora.care



Curtis Farringdon Account Manager 07939 842420 curtis.farringdon@accora.care



Nathan Smith Technical Sales Manager 01223 206100 nathan.smith@accora.care



Michael Sadler Product Specialist South West 07399 584409 michael.sadler@accora.care



Scott Wheeler Product Specialist Team Manager 07867 473198 scott.wheeler@accora.care



**Chris Seal** Product Specialist London & South East 07376 302557 chris.seal@accora.care



Mark Hodges Product Specialist Midlands 07584 560495 mark.hodges@accora.care



Mark Shires Product Specialist North 07984 606491 mark.shires@accora.care



Declan Moore Product Specialist Scotland 07584 071403 declan.moore@accora.care



Tino Paphitis Product Specialist London, Surrey, Kent, Essex 07376 302556 tino.paphitis@accora.care



## **Our services**

### **Joint visits**

We offer joint visits at your service user's home for the Accora range of products. Our Product Specialists recommend the minimum specification to meet the clinical need of your service user and work up from that as directed by the prescriber. If we believe that our products do not meet the clinical need, we'll advise the prescriber and suggest an alternative supplier.

### Sale or return

Accora products are available on a sale-or-return basis to enable users to trial a product before it is purchased. Sale-or-return provides complete peace of mind because it allows up to 14 days to evaluate accurately whether the product will meet the user's needs.

### **Customer service**

Accora's customer service team pride themselves on delivering an excellent customer service experience. Whether you need to request a quote, get an update on an order, enquire about a product or book a joint visit, they'll be pleased to help.

### **Professional support**

Since qualifying as an Occupational Therapist in 2006, Debra has worked with users with health and social care needs in community and inpatient settings.

In 2014 Debra split her time between her own case load of clients, her role with Accora and as equipment lead in the London Borough of Barnet equipment service, a role she continued until 2019.



Debra continues to practise as an occupational therapist alongside her role at Accora, which involves using her clinical knowledge and experience to inform product development, and to provide training from a clinical perspective on our product range. Debra uses her vast range of experiences as an Occupational Therapist and as equipment lead to support customers with client specific queries, and to help develop issuing criteria and sign off guidance for specialist equipment.

Contact Debra on 07967 702251 or debra.dunitz@accora.care



# **Accora Seating**

## Specialist seating and riser recliner chairs for community equipment stores...

- Easy adjustment of size, pressure and posture requirements before and after delivery
- Easy to recycle and reissue
- Reduces costs



### Overview

	Configura <sup>®</sup> Lite 2	Configura <sup>®</sup> Comfort	Configura <sup>®</sup> Rise & Recline
Tilt-in-space	<ul> <li>Image: A second s</li></ul>	✓	<ul> <li>Image: A second s</li></ul>
True tilt-in-space			
Riser function	<ul> <li>✓</li> </ul>	✓	<ul> <li>Image: A start of the start of</li></ul>
Flat rise		<i>✓</i>	🗸 (XL and S Tall)
Independent leg rest and tilt-in-space			
Backrest recline		<i>✓</i>	🗸 (in dual motor chair only)
Footplate			
Portable			
Depth adjustment options	✓	✓	<ul> <li>Image: A set of the set of the</li></ul>
Height adjustment options	<ul> <li>✓</li> </ul>	$\checkmark$	✓
Adjustable leg rest to accommodate height	<ul> <li>Image: A second s</li></ul>	✓	✓
Width adjustment options		<i>√</i>	
Pressure cushion options	<ul> <li>Image: A set of the set of the</li></ul>	✓	<ul> <li>Image: A second s</li></ul>
Posture management options	<ul> <li>Image: A second s</li></ul>	<i>√</i>	<ul> <li>Image: A set of the set of the</li></ul>
Transfer aid compatible		$\checkmark$	

Configura <sup>®</sup> Care	Configura <sup>®</sup> Advance	Configura <sup>®</sup> Porter	Configura <sup>®</sup> Bariatric
1	$\checkmark$	✓	✓ (Tilt-in-space version)
1	$\checkmark$	$\checkmark$	1
		✓	1
			1
1	$\checkmark$	✓	✔ (Tilt-in-space version)
	✓	1	1
<ul> <li>✓</li> </ul>	$\checkmark$	✓	
1	✓	✓	
<ul> <li>✓</li> </ul>	$\checkmark$	✓	✓
	$\checkmark$	<ul> <li>(16" and 18" versions only)</li> </ul>	✓ (16" and 18" versions only)
	$\checkmark$		
1	✓	<i>s</i>	1
<ul> <li>✓</li> </ul>	✓	1	1
1	✓	<i>√</i>	✓
	✓		

### Not sure what chair to use?

Call our team on 01223 206100 or email info@accora.care to book a joint visit.


### Chair accessories

			Configura <sup>®</sup> Lite 2	Configura <sup>®</sup> Comfort
	1" Height Adjustment Kit	<b>.</b>		1
Size	2" Height Adjustment Kit	-	$\checkmark$	1
	Depth Adjustment Kit	AT I	$\checkmark$	1
	Removable Cushion Core		$\checkmark$	1
	Castellated Cushion option			
Pressure	Cushionair option		$\checkmark$	<ul> <li>Image: A second s</li></ul>
	Drop-in Overlay options			
	Visco Cushion		$\checkmark$	1
	Waterfall Backrest (std)	Z	$\checkmark$	<ul> <li>Image: A second s</li></ul>
	Lateral Support Wedges		$\checkmark$	1
	Lateral Support Backrest			1
	Adjustable Lateral Support			<ul> <li>Image: A start of the start of</li></ul>
	Profiled Headrest			<ul> <li>Image: A second s</li></ul>
	Drop-down Armrest	The second		1

Configura <sup>®</sup> Rise & Recline	Configura <sup>®</sup> Care	Configura <sup>®</sup> Advance	Configura <sup>®</sup> Porter	Configura <sup>®</sup> Bariatric
1				
✓				
✓	1	1	1	1
	1	<i>✓</i>		
✓	1			
✓	1	1	1	1
1			1	1
✓	1	1	1	1
✓	1	1	1	1
✓	1	1	1	1
✓	1	1	1	
✓	$\checkmark$	1	<b>√</b>	
✓	$\checkmark$	1	$\checkmark$	1
		<i>✓</i>		

### Configura® Concept

The Configura chair range by Accora introduces a revolutionary way of issuing seating which is saving many departments time, effort and money. Here's how it works:



The new Configura Lite 2 is Accora's latest riser recliner which combines many of the benefits of the Configura Comfort chair, together with the cost-effectiveness of the previous Configura Lite.



Chair size	e Part number	Seat depth	Seat height	Seat width	Weight capacity
Small	CHAIR-0-CM2-140	450 mm (18")	400 mm (16")	450 mm (18")	25 stone
Sinui	SITION CHAIR-U-CM2-140	500 mm (20")	450 mm (18")**		
Medium	CHAIR-0-CM2-240	450 mm (18")	450 mm (18")	500 mm (20")	25 stone
Medium	CHAIR U-CIVI2-240	500 mm (20")	500 mm (20")**	500 min (20 )	20 StOFIE

\*\* 2" Height Adjustment Kit required

Standard size

### **Configura**<sup>®</sup> Comfort

### Configura<sup>®</sup> Comfort features and accessories

### Configura Comfort is a revolutionary chair providing an unprecedented level of functionality and cost-effectiveness.

'Reduced shear' recline to maximise pelvic stability

Integrated self-supply overlay for easy integration of existing cushions.

Compatible with transfer aids including Cricket, ReTurn 7500 and Molift Raiser

Braked castors

Chair size

Small

Medium

Larae



375 mm (15")°

400 mm (16")

425 mm (17")\*

450 mm (18")\* 425 mm (17")°

450 mm (18")

475 mm (19")\*

500 mm (20")\* 425 mm (17")

450 mm (18")

475 mm (19")\*

500 mm (20")\*

450 mm (18")

500 mm (20")

550 mm (22")

450 mm (18")

500 mm (20")

450 mm (18")

500 mm (20")

450 mm (18")

500 mm (20")

Pressure reducing fabric on backrest and armrest as standard

> Colour-matched vapour permeable fabric

Knock-down design for easier storage and delivery

Vinyl fabric option

25 stone

25 stone

25 stone

01223 206100





Braked castors

Integrated self-supply overlay. Aperture dimensions 17.5" x 17.5" x 3"











Pressure-reducing fabric

on backrest and armrest

as standard





Compatible with transfer

aids including Cricket, ReTurn

7500 and Molift Raiser

Colour-matched vapour-permeable fabric



'Reduced shear' recline to

maximise pelvic stability

Drop-down armrest option




More depth adjustment

options using backrest

brackets. Reduces the need

for depth adjustment kits



Vinyl fabric option

\*\* Includes 2" Height Adjustment Kit

- Includes 1" Height Adjustment Kit
- The castors are removed to make the chair lower

CHAIR-0-CM1-140

CHAIR-0-CM1-240

CHAIR-0-CMI-340

Standard size

Seat depth is adjustable by 50 mm  $(2^{\circ})$  as standard without extra parts

### **Configura**<sup>®</sup> Bariatric

Configura Bariatric chair offers unprecedented adjustability, flexibility and peace of mind for the prescriber. There are three chairs in this range: 40 stone dual-motor chair, 40 stone dual-motor tilt-in-space chair, and 50 stone dual-motor chair.



450 mm (18") 700 mm (28") 400 mm (16") CFB282018DT 500 mm (20") 650 mm (26") 40 stone 40 stone dual-motor 450 mm (18") 550 mm (22") 600 mm (24") 450 mm (18") 400 mm (16") 700 mm (28") 40 stone dualmotor tilt-in-space CFBTIS282018DT 500 mm (20") 450 mm (18") 650 mm (26") 40 stone (with independent legrest) 550 mm (22") 500 mm (20") 600 mm (24") 450 mm (18") 850 mm (34") CFBSUP303234DT 500 mm (20") 450 mm (18") 800 mm (32") 50 stone 50 stone dual-motor 550 mm (22") 750 mm (30")

° Depth Adjustment kit (2" decrease in seat depth)

<sup>∞</sup> 2 x Depth Adjustment Kits (4" decrease in seat depth)

Standard size

Configura Porter is a revolutionary product which allows a level of functionality not previously possible with specialist seating. For clients with fluctuating functional levels.



Chair size	Part number	Seat depth	Seat height	Seat width	Weight capacity
		450 mm (18") <sup>~~</sup>		400 mm (16")	
Small	CFP182016DT	500 mm (20") °	400 mm (16")	450 mm (18")	20 stone
		550 mm (22")		500 mm (20")	
		450 mm (18") <sup>∞</sup>		400 mm (16")	
Medium	CFP182020DT	500 mm (20") °	450 mm (18")	450 mm (18")	20 stone
		550 mm (22")		500 mm (20")	

• Depth Adjustment kit (2" decrease in seat depth)

<sup>oo</sup> 2 x Depth Adjustment Kits (4" decrease in seat depth)

Standard size

### Configura® Advance

### **Configura**<sup>®</sup> Advance features and accessories

A comprehensive posture management system designed specifically for high need users.



Part number	Adjustment	Range	Weight capacity
	Height	330mm (13") – 530mm (21")	
CHAIR-0-SC1-030	Depth	400mm (16") – 550mm (22")	25 stone
	Width	350mm (14") – 550mm (22")	









Integrated self-supply overlay, aperture dimensions 17.5" x 17.5" x 3"

Continuous overlay

for thigh support



Full-length compensation on legrest, giving superior leg elevation

Legrest elevation with 37

degrees of tilt-in-space



Angle height-adjustable footplate



Forward tilt (option on all chairs)



Dropdown armrest – helps with transfers and sling insertion

### Configura® Rise & Recline

### Configura® Care

With six standard sizes and 36 configurations, Rise & Recline is one of the most versatile chairs in the Configura range. It's appropriate for ambulant, semi-ambulant and non-ambulant users, accommodation and encouragement of good posture, and users at risk of pressure damage.

High-risk pressure reduction as standard

Wide range of size requirements

36 size options

Lead time 2-3 days



8 posture options

Single, dual-motor tilt-in-space or dual-motor non-tilt-in-space

4 integrated, interchangeable pressure care options

Chair size	Part number	Seat depth	Seat height	Seat width	Weight capacity									
		400 mm (16")°	375 mm (15")											
Extra Small	CFXSD161815DT	( )	400 mm (16")*	400 mm (16")	20 stone									
		450 mm (18")	425 mm (17")**											
		450 mm (18")°	450 mm (18°)											
Small Tall	CFSTD182018DT		475 mm (19")*	450 mm (18")	25 stone									
		500 mm (20")	500 mm (20")**											
	Medium Low CFMLP202016DT	450 mm (18")°	400 mm (16")											
Medium Low			425 mm (17")*	500 mm (20")	20 stone									
												500 mm (20")	450 mm (18")**	
	CFLLP222016DT	CFLLP222016DT	CFLLP222016DT	450 mm (18")° 400 mm (16")										
Large Low				CFLLP222016DT	CFLLP222016DT	CFLLP222016DT	CFLLP222016DT	CFLLP222016DT	CFLLP222016DT	CFLLP222016DT		425 mm (17")*	550 mm (22")	20 stone
		500 mm (20")	450 mm (18")**											
		450 mm (18")°	450 mm (18°)											
Extra Large	CFXLD242018DT	475 mm (19")* 600 mm (24")		600 mm (24")	25 stone									
		500 mm (20")	500 mm (20")**											
			400 mm (16")											
Extra Large Low CFXLLP24	CFXLLP242016DT		425 mm (17")*	600 mm (24")	20 stone									
		500 mm (20")	450 mm (18")**											

\* 1" Height Adjustment Kit required 👋 2" Height Adjustment Kit required 🛸 2" Depth Adjustment Kit required 📕 Standard size

Configura Care is an affordable and adjustable chair. The Configura Care chair is easy to recycle, adjustable to suit the user's changing needs and cost-effective – giving better budget management.



Chair size	Part number	Seat depth	Seat height	Seat width	Weight capacity	
	400 mm (16") <sup>∞</sup>		300 mm (12")			
Small	Small CFC240510DT	450 mm (18")°	300 – 500 mm	350 mm (14")	20 stone	
Sindi		( - )	000 00011111	400 mm (16")	20 3:0116	
		500 mm (20")		450 mm (18")		
	Medium CEC450510DT		450 mm (18") <sup>∞</sup>		400 mm (16")	
Medium		500 mm (20")°	300 – 530 mm	450 mm (18")	25 stone	
Medium CFC450510D1	500 mm (20 )	300 33011111	500 mm (20")	zo stone		
		550 mm (22")		550 mm (22")		

° Includes 2" Depth Adjustment Kit

<sup>oo</sup> Includes 2 x 2" Depth Adjustment Kits

Standard size

Please note: seat-to-floor height is 560 mm on Medium and 530 mm on Small. The seat-to-footplate height can be adjusted from 300 mm to 530 mm on Medium and 400 mm to 500 mm on Small. The footplate is an option.

### Configura® Custom

Configura Custom chair combines the functionality of Configura Comfort, with the added benefit of a range of fabrics to choose from.





Chair size	Part number	Seat depth	Seat height	Seat width	Weight capacity
		450 mm (18") 375 mm (15")	375 mm (15")°		
Small	CHAIR-0-CM1-140	450 11111 (16.)	400 mm (16")	450 mm (18")	25 stone
Smail	CHAIR O CIVILIAO	500 mm (20")	425 mm (17")*	430 11111 (10 )	20 310116
		500 mm (20 )	450 mm (18")**		
	Medium CHAIR-0-CMI-240	450 mm (18")	425 mm (17")°	500 mm (20")	25 stone
Modium			450 mm (18")		
Mediditti		500 mm (20")	475 mm (19")*		
		500 (1111 (20 )	500 mm (20")**		
		450 mm (18")	425 mm (17")°		
Large CHAIR-0-CM1-340	400 mm (10 )	450 mm (18")	550 mm (22")	25 stone	
	500 mm (20")	475 mm (19")*	000 mm (22)	20 310116	
		500 mm (20") 500	500 mm (20")**		

" Includes 2" Height Adjustment Kit

Includes 1" Height Adjustment Kit

° The castors are removed to make the chair lower

Standard size

Seat depth is adjustable by 50 mm (2") as standard without extra parts



## **Accora Beds**

#### www.accora.care

## Bed fall, or bed roll?

- Reduce the risk of injury in the event of a user fall
- Reduce the risk of care giver injury
- Enable independence and mobility






### **FloorBed**<sup>®</sup> ]

FloorBed 1 has a mattress platform height of just 7.1 cm, greatly reducing the risk of injury from falls, and can help to promote independence. It rises to nursing height to facilitate care tasks and offers comprehensive functional and clinical benefits to maximise user safety.



Overall dimensions	925 x 2255 mm (36.4" x 88.8")			
Mattress size	900 x 2000 mm (35.4" x 78.7")			
Safe working load	185 kg (29 st, 408 lbs)			
Patient weight	150 kg (24 st, 330 lbs)			
Lifting range	71 – 650 mm (2.8" – 25.6") (excl. mattress)			
Profiling functions	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg			
Min. height of service user	146 cm (57.5°)			
Part number	NSB-0-FLI-200			

### FloorBed 2 lowers to just 7.1 cm and rises to 80 cm, enabling carers to complete care tasks at an appropriate height.



Overall dimensions	935 x 2285 mm (36.8" x 90.0")
Mattress size	900 x 2000 mm (35.4* x 78.7*)
Safe working load	185 kg (29 st, 408 lbs)
Patient weight	150 kg (24 st, 330 lbs)
Lifting range	71 - 800 mm (2.8"- 31.5") (excl. mattress)
Profiling functions	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg
Min. height of service user	146 cm (57.5")
Part number	NSB-0-FL2-200

### **Floorbed**<sup>®</sup> accessories

### All FloorBeds are available with a range of accessories, detailed below:





Integrated padded side rails



Safety sleeve reducing risk of entrapment as the platform mechanism is protected





Bed lever to assist with bed mobility



High safety mat, 22 cm high.

Same height as FloorBed

with 15 cm foam mattress

Rotating bed lever for extra support in movement



Lifting pole to assist with bed mobility



Mattress infill



Bed extension - increasing mattress platform length to 2200 mm

## **Junior Beds**

- Appropriate for service users from 75 cm tall
- Meets the BS EN 50637 standard, helping to keep your service user safe, and you compliant









### FloorBed® 1 Junior

Accora have developed a simple, safe and cost-effective solution to enable departments to stay compliant with the new BS EN 50637 standard for children's beds. FloorBed 1 Junior is suitable for service users with a height of 75 cm and upwards.



Overall dimensions	925 x 2255 mm (36.4" x 88.8")
Mattress size	900 x 2000 mm (35.4" x 78.7")
Safe working load	185 kg (29 st, 408 lbs)
Patient weight	150 kg (24 st, 330 lbs)
Lifting range	71 – 650 mm (2.8" – 25.6") (excl. mattress)
Profiling functions	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg
Min. height of service user	75 cm (29.5")
Part number	NSB-0-FLI-200 + JUNKIT-0-FLI-000

#### Features:



Built-in and lockable pocket for storing the handset securely to prevent tampering

### **Accessories:**



FirmEdge mattress – designed with a firm edge to further reduce risks with entrapment and allow easier transfers. The firm edge helps service users orientate themselves in bed, and can make transfers easier



Head and foot board

bumpers reduce impact

and minimise the risk

of injury

Side rails with integrated

bolster made of soft

materials



Safety sleeve reducing risk of entrapment as the platform mechanism is protected



Detailed image of the side rail with built-in side bolster


### FloorBed 1 Junior

Comprises	Part number	Descr	iption
FloorBed 1	NSB-0-FL1-200	FloorBed 1	
Junior Kit	JUNKIT-0-FLI-000	Comprises Head and Foot Board Bumpers (BMHNE-0-FL4-000) and Safety Sleeve (SLE-0-FL4-000)	

PLEASE NOTE: FloorBed 1 Junior can only meet the new children's standard when the installation of the full **Junior Kit** is complete, as specified in the table.

### FloorBed 1 Junior accessories

Comprises	Part number	Descr	iption
Mattress	VISMATB-0-FM2-000	FirmEdge <b>Allevia</b> ® Green Mattress	11 e
Side Rail with Side Bolster	SDR-0-FL4-000	Side Rails for <b>FloorBed</b> 1 Junior	



### Quick guide Children's standard



### What is it?

The new children's standard (BS EN 50637:2017 – Medical electrical equipment – Particular requirements for the basic safety and essential performance of medical beds for children) is a new standard for medical beds suitable for children.

### Why has it been introduced?

Adult beds have an intended service user group which is defined as below (adult standard is BS EN 60601-2-52):

- Over 40 kg
- Over 1.46 m in height
- With a BMI of more than 17

Anyone outside these dimensions, particularly children, therefore fall out of this user group.

### What is different in the new standard?

The main difference in the standard is the size of the gaps between side rails, mattress platform and closing gaps. These are much smaller than with adult beds. There are also extra tests that must be conducted on the bed to enable it to meet the standard, that are not required on the adult bed.

### How will it affect me?

Manufacturers must meet these standards to be able to state that their beds are suitable for children. For Local Authorities, this means that when issuing beds for children, or those who are outside of the intended service user groups for adults, an appropriate bed is used.

### What has Accora done to meet the standard?

Accora have developed a safe, simple, costeffective solution for our beds to enable community equipment stores / prescribers to be able to comply with the new BS EN 50637 standard. Please contact us for more information and we'll be happy to discuss this in further detail with you.

## A low bed that goes to a full nursing height

- Lifts as high as a conventional profiling bed (80 cm)
- Goes as low as a conventional low bed (22 cm)
- Compliant with BS EN 60601-2-52



### CommunityBed

CommunityBed allows equipment services to combine two beds into one. By combining the lifting height of a conventional profiling bed (800 mm) and the lowness of a conventional low bed (220 mm), CommunityBed allows equipment services to simplify their range and reduce stock. CommunityBed is fully compliant with BS EN 60601-2-52 and compatible with a broad range of accessories including integrated, padded-mesh side rails.



Overall dimensions	930 x 2250 mm (36 x 89")
Mattress size	900 x 2000 mm (36 x 79*)
Safe working load	215 kg (34 st)
Patient weight	180 kg (28 st)
Lifting range	220 – 800 mm (8.6" – 31.5")
Profiling functions	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg
Min. height of service user	146 cm (57.5")
Part number	NSB-0-CB1-000

#### **Features:**





Easy to use handset

Fully profiling



Excellent height for carers



Lowest position at 22 cm

### **Accessories:**









itegrated side rails <sup>.</sup> Iowered



## **Pressed for time?**

### **Book a free Accora product awareness session**

These sessions will raise awareness of our **FloorBed** and **Configura** ranges. Enjoy a free lunch, courtesy of Accora, at the same time.

Typical learning outcomes from a lunch and learn session include:

- How to issue appropriate seating
- How to reduce risk of bed fall injuries
- Understanding new children's bed regulations

Certificate of attendance issued. To book a session, please contact us using the details below.

## Notes

### Notes

Notes	Notes

Notes	Notes

	-
Notes	
	_
	-
	-
	-
	_
	_
	-
	-
	_
	_
	_
	_
	_
	-
	-
	-
	_

01223 206100

Charter House, Barrington Road Orwell, Cambridge SG8 5QP

- **T:**+44 (0)1223 206100
- E: info@accora.care
- W: www.accora.care

and ™ are trademarks of Accora Limited.
 As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.